


XL3000i Autoclavable Pipettes Spare Parts List 2012



	Catalog Number	Description	Pipette	Price
1	SP19360 SP19361 SP19362 SP19363 SP19364 SP19365 SP19366 SP19367	Pushbutton Pushbutton Pushbutton Pushbutton Pushbutton Pushbutton Pushbutton Pushbutton	XL2 XL10 XL20 XL100 XL200 XL1000 XL5000 XL10000	\$9.36
1A	SP19349	Knob (lower Portion of Pushbutton)	XL	\$5.46
2	SP001-B	Upper arbor	XL	\$32.76
3	SP003LM	Calibration knob Complete	XL	\$42.00
3a	SP003-C	Calibration screw	XL	\$6.86
5	SP003-E	Fixing screw	XL	\$1.56
6a	SP9385	O-ring	XL	\$1.87
7	SP009-XL	Shaft nut	XL2-1000	\$3.12
7*	SP009-XL2	Shaft nut	XL2-1000	\$7.80
8	SP005-2 SP005-10 SP005-20 SP005-100 SP005-200 SP005-1000 SP005-5000 SP19463	Piston assembly Piston assembly Piston assembly Piston assembly Piston assembly Piston assembly Piston assembly Piston assembly	XL2 XL10 XL20 XL100 XL200 XL1000 XL5000 XL10000	\$40.56 \$40.56 \$40.56 \$40.56 \$40.56 \$40.56 \$40.56 \$65.00
8*	SP19624 SP19398 SP19399 SP19400 SP29301 SP29302 SP9275 SP19482	Piston Assembly Piston Assembly Piston Assembly Piston Assembly Piston Assembly Piston Assembly Piston Assembly Piston Assembly	XL2 XL10 XL20 XL100 XL200 XL1000 XL5000 XL10000	\$40.56 \$40.56 \$40.56 \$40.56 \$40.56 \$40.56 \$40.56 \$60.00
9	SP005A-100	Spring	XL2, XL10, XL20, XL100	\$3.12
9*	SP29303	Spring	XL2, XL10, XL20, XL100	\$3.12
10	SP005A-200	Spring	XL200	\$3.64
11	SP005D	Stopper bushing	XL2, XL10, XL20 XL100, XL200	\$2.00

	Catalog Number	Description	Pipette	Price
12	SP9321	Compression spring	XL2, XL10, XL20 XL100, XL200	\$3.12
13	SP005B	Bush	XL2, XL10, XL20	\$2.00
13*	SP29306	Bush	XL100, XL200	\$4.68
14	SP011-2 SP011-10 SP011-20 SP011-100 SP011-200 SP011-1000 SP011-5000 SP19464	Compression bushing Compression bushing Compression bushing Compression bushing Compression bushing Compression bushing Compression bushing Compression bushing	XL2 XL10 XL20 XL100 XL200 XL1000 XL5000 XL10000	\$4.68 \$3.64 \$3.64 \$3.64 \$3.64 \$3.64 \$3.64 \$3.64
15	SP006-2 SP006-10 SP006-20 SP006-100 SP006-200 * SP006-200A* SP006-1000 SP006-5000 SP19465	Teflon seal Teflon seal Teflon seal Teflon seal Teflon seal Teflon seal Teflon seal Teflon seal Teflon seal	XL2 XL10 XL20 XL100 XL200 XL200 XL1000 XL5000 XL10000	\$6.24 \$6.24 \$6.24 \$6.24 \$6.24 \$6.24 \$6.24 \$6.24 \$6.24
16	SP007-12 SP007-10 SP007-100 SP007-200 * SP007-200A* SP007-1000 SP007-5000 SP19002	O-ring 0.74x1.02 O-ring O-ring O-ring O-ring O-ring O-ring O-ring O-ring 25x2	XL2 XL10, XL20 XL100 XL200 XL200 XL1000 XL5000 XL10000	\$9.00 \$4.16 \$4.16 \$4.16 \$4.16 \$4.16 \$4.16 \$4.16
17	SP010-10XL SP9442 SP010-100XL SP010-200XL SP010-1000XL SP19452 SP19453	Tip ejector Tip ejector Tip ejector Tip ejector Tip ejector Tip ejector Tip ejector	XL2, XL10 XL20 XL100 XL200 XL1000 XL5000 XL10000	\$19.76 \$19.76 \$19.76 \$19.76 \$19.76 \$31.20 \$42.00
17*	SP19446 SP19447 Sp19448 SP19449 SP19450 SP29418	Tip ejector Tip ejector Tip ejector Tip ejector Tip ejector Tip ejector	XL2, XL10 XL20 XL100 XL200 XL1000 XL1000 **	\$18.00 \$18.00 \$18.00 \$18.00 \$18.00 \$18.00

	Catalog Number	Description	Pipette	Price
18	SP008-2LM	Shaft	XL2	\$14.56
	SP008-10LM	Shaft	XL10	\$14.56
	SP008-20LM	Shaft	XL20	\$14.56
	SP008-100LM	Shaft	XL100	\$14.56
	SP008-200LM	Shaft	XL200	\$14.56
	SP008-1000LM	Shaft	XL1000	\$14.56
	SP008-1000LM2	Shaft	XL1000**	\$14.56
	SP008-5000LM	Shaft	XL5000	\$14.56
	SP008-10000LM	Shaft	XL10000	\$14.56
20	SP003-A	X-ring	XL	\$3.12
21	SP012	Pipette calibration service	XL	\$50.00
22	SP002-EB1	Ejector button I	XL	\$6.24
22a	SP002-EB2	Ejector button II	XL	\$6.24
25	SP9438	Ejector arbor	XL2-1000	\$18.72
	SP19402	Ejector arbor	XL5000, XL10000	\$18.72
25a	SP9397	Ring 3z	XL	\$2.10
26	SM9479	Calibration key	XL	\$12.48
27	SP19454	Filter	XL5000, XL10000	\$1.04
28	SP19466	Ejector cap	XL5000	\$9.60
	SP19469		XL10000	

* Item for new (low forces) version only